

Overview

SINGLE-USE HYSTEROSCOPE

MacroLux Seleneview and ViewHub

Order Code	Product description	Included in delivery
H035 	Single Use Operative Hysteroscope Outer diameter: 3.5 mm Working channel diameter: 1.8 mm Resolution: 160k Pixel	Box Quantity: 10 Units, delivered sterile
H042 	Single Use Operative Hysteroscope Outer diameter: 4.2 mm Working channel diameter: 2.0 mm Resolution: 920k Pixel	Box Quantity: 10 Units, delivered sterile
HC29 	Single Use Diagnostic Hysteroscope Outer diameter: 2.9 mm Resolution: 160k Pixel	Box Quantity: 10 Units, delivered sterile
HC35 	Single Use Diagnostic Hysteroscope Outer diameter: 3.5 mm Resolution: 920k Pixel	Box Quantity: 10 Units, delivered sterile
S13 	Video Processor 13.3 inches Including Memory Card, Cardreader and Power Supply Type: HD touch screen Resolution: 1920 x 1080 pixels	Box Quantity: 10 Units, delivered sterile