






## ATOMS System

Order Code		Product	Included in delivery
<b>ATS5041</b>		<p><b>ATOMS</b> Adjustable system for treatment of male stress urinary incontinence System contains:</p> <ul style="list-style-type: none"> <li>• implant with soft, adjustable suburethral cushion</li> <li>• integrated, silicone-coated titanium port</li> <li>• non-coring port needles</li> </ul>	<p>Single-use system, delivered sterile</p> <p>For the ATOMS surgery with ATS5041 the following product is needed: 1x TOA5130 TOA Tunneller or 1x ATS5010 ATOMS Tunneller</p>
<b>TOA5130</b>		<p><b>A.M.I. TOA Tunneller</b> Reusable instrument of stainless steel for outside-in tunnelling. Treatment of male stress urinary incontinence with the A.M.I. ATOMS System</p>	<p>1 pair, helical (left and right), delivered non-sterile steam autoclavable</p>
<b>ATS5010</b>		<p><b>ATOMS Tunneller</b> Reusable instrument of stainless steel for outside-in tunnelling. Treatment of male stress urinary incontinence with the A.M.I. ATOMS System</p>	<p>1 pair, helical (left and right), delivered non-sterile steam autoclavable</p>
<b>ATS5051</b>		<p><b>Scrotal Port for ATOMS</b> Scrotal port with catheter, tubing connector and 2 port needles</p>	<p>1 unit, delivered sterile</p>
<b>ATS5061</b>		<p><b>Tubing Connector</b> Connection piece made from titanium for catheter-to-catheter connection, e.g. for port-revision with ATOMS implant. Includes catheter spare piece of 200 mm</p>	<p>2 units in 1 box, delivered sterile</p>
<b>ATS5031</b>		<p><b>Tubing Plug</b> Made from titanium For temporary closure of ATOMS catheter</p>	<p>2 units in 1 box, delivered sterile</p>
<b>SFN 0930 G</b>		<p><b>Port Needle</b> Port needle for adjustment for the ATOMS Scrotal Port Needle Length: 30 mm, 20 G x 30 mm</p>	<p>25 needles / box, delivered sterile</p>