












Overview


Treatment of Hemorrhoids

Order Code		Product description	Included in delivery
TRI2010		<p>Wi-3 HAL-RAR Unit System for HAL-RAR surgery with integrated bluetooth technology.</p> <p>1 system consists of:</p> <ul style="list-style-type: none"> • Doppler-ultrasound electronic unit with integrated bluetooth submitter • TRI2020 Wi-3 Battery Cap • TRI2030 Wi-3 Fixation nut • TRI2040 Wi-3 Battery Set 	Delivered non-sterile, steam autoclavable Multi-use (50 applications) Batteries: 2 units, AA NiMH, 2100 mAh
TRI2110		<p>Speaker for TRILOGY Bluetooth speaker for the optimal reproduction of the TRILOGY Unit's signal (TRI2010) during HAL-RAR surgery</p>	1 unit, supplied with USB cable
TRI2210		<p>Adaptor for Speaker USB adaptor for connecting the speaker Manufacturer Frimo FOX12-X11-USB</p>	1 unit, delivered non-sterile
TRI2220		<p>Connector Europe Replaceable connector for connection to adaptor, Manufacturer Frimo</p>	1 unit, delivered non-sterile
TRI2020		<p>Wi-3 Battery Cap Reusable battery cap for the TRILOGY Unit (TRI2010)</p>	1 unit, delivered non-sterile, steam autoclavable
TRI2030		<p>Wi-3 Fixation Nut Reusable fixation nut for the fixation of probes onto the TRILOGY Unit (TRI2010)</p>	1 unit, delivered non-sterile, steam autoclavable
TRI2040		<p>Wi-3 Battery Set 1 package consisting of 2 rechargeable batteries (AA, NiMH, 2100 mAh)</p>	2 pieces, delivered non-sterile
TRI2070		<p>Wi-3 Arm Reusable holding wheel to enhance holding position for the assistance during surgery</p>	1 unit, delivered non-sterile, steam autoclavable
RAR2081		<p>RAR Flexi Probe Disposable probe and sleeve set for performing HAL and RAR procedures. Probe with asymmetric design for the gradual release of mucosa</p>	1 box, delivered sterile 5 x RAR Flexi probe and sleeve
RAR2181		<p>HAL-RAR Procedure Kit – Flexi</p>	1 kit, delivered sterile 1 x RAR Flexi probe and sleeve 8 x A.M.I. HAL Suture 1 x Single-use knotpusher designed especially to fit the A.M.I. probes
AHAL 70		<p>A.M.I. HAL Suture Suture material for HAL and RAR procedures.</p>	1 box, delivered sterile 36 x sutures, 5/8 needle absorbable, 2/0, 75cm

<p>AHN 006</p>		<p>A.M.I. HAL Needleholder Stainless steel needleholder designed especially to fit the A.M.I. probes</p>	<p>1 instrument, delivered non-sterile, steam autoclavable</p>
<p>AHK 007</p>		<p>A.M.I. HAL Knotpusher Stainless steel knotpusher to facilitate knot-tying inside the A.M.I. probes</p>	<p>1 instrument, delivered non-sterile, steam autoclavable</p>




Overview

Fecal Incontinence

Order Code	Product description	Included in delivery
<p>AAB 300</p> 	<p>Soft Anal Band System for treatment of fecal incontinence.</p>	<p>1 system, delivered sterile 1 x Soft Anal Band Ring M/L 1 x Soft Anal Band Activator M/L 1 x Soft Anal Band Valve 1 x Soft Anal Band Port 1 x Soft Anal Band Port Connector 1 x Soft Anal Band Ring Connector</p>

Overview

Anal Abscesses & Fistulas

Order Code	Product description	Included in delivery
<p>CDS2001</p> 	<p>Comfort Drain System consisting of a flexible stainless steel probe for the location of fistula passages and a unique knotless 360° lockable silicone drain with a polypropylene connector. System suitable for long-term drainage of anal fistulas.</p>	<p>1 box, delivered sterile 10 x set</p>
<p>CDS2010</p> 	<p>Comfort Drain Grasper Reusable instrument to facilitate the closure of the Comfort Drain CDS2001.</p>	<p>1 box, non-sterile, steam autoclavable 1 x Comfort Drain Grasper</p>
<p>FIX2001</p> 	<p>FiXcision Enables the circumferential tissue removal to treat simple anal fistulas.</p> <p>Set consisting of probe to discover the fistula tract, guide for tissue compression around the fistula tract, base plate and circular blade for precise excision of the fistula tract.</p>	<p>1 box, delivered sterile 10 x set</p>